



/ CELL-DISRUPTION / DEAGGLOMERATION / DISPERSIONS

/ EMULSIONS / LIPOSOMES / PARTICLE SIZE REDUCTION

/ POLYMER MW REDUCTION / STABILIZATION / ...AND MORE



## **PSI-300**

- up to 2068 bar
- / up to 50Hz 500 L/h 60Hz 550 L/h
- virtually unlimited throughput in modular configuration
- / less than 70 dBA
- / flexible configurations

Thanks to its flexible architecture and efficient hydraulic system, the **PSI-300** can achieve best-in-class throughput in a compact space. It can be customized to meet the most demanding requirements, and optimized to the workload to provide the most cost-effective solution.

## **FOOTPRINT**

W×L×H 1310×1900×2245 mm

## **POWER**

**50Hz**3×400V 50Hz

3×480V 60Hz

P<sub>N</sub> 25 kW, I<sub>N</sub> 56A

P<sub>N</sub> 30 kW, I<sub>N</sub> 56A

**WEIGHT** 1850 kg

**FLOWRATE** 

**50Hz** 500 L/h **60Hz** 550 L/h

**MAX PRESSURE** 

2068 bar

The choice of the best materials and finishing, the long list of options, make the PSI-300 suitable for every industry, including the most demanding pharmaceutical, biochemical and food processing environments.

Several configurations available:

up to 50Hz 300 L/h 60Hz 330 L/h @ 2068 bar up to 50Hz 400 L/h 60Hz 440 L/h @ 1500 bar up to 50Hz 500 L/h 60Hz 550 L/h @ 1200 bar

Optional CIP and SIP, dedicated enclosures for specific products, ATEX option, additional process sensors, data acquisition and recording.

## **CONTACT US**

info@psi-instruments.com www.psi-instruments.com RIFFRIMENTO PER L'ITALLA



Qi srl t +39 06 9105461 www.qitech.it | SalesQi@qitech.it

